



Simplicity+		inch	mm
A	Minimum swivel radius	26½	655
B	Minimum wall to seat back	3	75
C	Seat depth	16¾	426
D	Minimum folded width to footplate	15¼	385
E	Minimum open width to edge of footplate	25½	650
F	Armrest width – external	21¾	550
G	Armrest width – internal	18	460
H	Seat back height	17¾	453
I	Footplate to seat height – standard seat option – perch seat option	18¼-21¼ 26¾-34½	465-540 677-877
J	Minimum footplate height	3¼	85
Bottom parking space required		21½ - 26½ 550mm (56°) min to 675 at shallow angle (27°)	
Minimum track intrusion into staircase		4½	115
Minimum track overhang at top of stairs		3½	150
Minimum track intrusion at bottom of stairs		15¾	250
Maximum weight limit – heavy duty version		302lb / 160.45 kgs 350lb / 159.09 kgs	
Speed		4¾"/sec (approx)	
Power supply		24v DC battery	
Charger		100-240v (50-60Hz)	
Maximum staircase angle		56° (52° with hinge)	



Sterling Stairlifts from Handicare are backed by a 2 year parts warranty plus a LIFETIME warranty on the gearbox and motor.

Sterling Stairlifts - Handicare Accessibility
 10 Industrial Park Drive
 Mount Pocono, PA 18344
 Toll Free: 866-276-5438
 Local: 570-895-2480
 Fax: 570-895-2481
info@sterling-stairlifts.com www.sterling-stairlifts.com



Please be aware that actual colors may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



A perch option is also available on the 1000



Minivator 1000
 Straight stairlift

Charron
 Medical Services
www.nhlift.com
 1-800-660-7221

The benefits of a Sterling 1000 from Charron Medical Services



If you are thinking about installing a stairlift the Sterling 1000 offers you a safe and reliable solution.

Handicare's Sterling Stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Sterling 1000 stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

With a choice of six upholstery colors the Sterling 1000 will compliment most color schemes.

The extremely slim track is fitted to the stairs, not the wall, so you will not need to

redecorate following the installation of your Sterling 1000 stairlift.

Like all Handicare stairlifts the Sterling 1000 has a battery backup which means it will continue to work in the event of a power outage so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the slim track maximises clear space on the stairway. It can be supplied with a manual or powered hinged track if required.

The Sterling 1000 also comes with two remote controls which means, once folded, you can easily send it to the upstairs charge point freeing up space in your hallway.

The Sterling 1000 is also available in a special outdoor model.



Safety As Standard

Your Sterling 1000 stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

The Sterling 1000 stairlift has an easy release retractable seat belt for added peace of mind.

Smooth & Quiet

Because of the high quality motors we use in our Minivator 1000 stairlifts we are able to offer perhaps the quietest and smoothest stairlift ride available. All backed by a lifetime gearbox and motor guarantee.

Easy to Use

The operating toggle is easy to use. Simply hold it towards the direction you wish to travel, up or down, and your Sterling 1000 will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.

